**3M**

**Petrifilm™**

6490/6491

---

-  (GB) **Staph Express Count System**
-  (FR) **Système de numération Staph Express**
-  (DE) **Staph Express Zählsystem**
-  (IT) **Sistema per il conteggio Staph Express**
-  (ES) **Sistema de recuento Staph Express**
-  (NL) **Staph Express Count System**
-  (SE) **Staph Express beräkningssystem**
-  (FI) **Staph Express Count System -pikakasvatusalusta**
-  (PT) **Sistema de Contagem Expressa de Staphylococcus**
-  (GR) **Σύστημα Ταχείας Καταμέτρησης Staph**
-  (JP) **黄色ブドウ球菌エクスプレス測定システム**

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#### AOAC Official Method 2003.07

3M Petrifilm Staph Express Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Types of Processed and Prepared Foods

#### AOAC Official Method 2003.08

3M Petrifilm Staph Express Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Dairy Foods

#### AOAC Official Method 2003.11

3M Petrifilm Staph Express Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Meat, Seafood, and Poultry

For specific details, refer to the specific AOAC® Official Method<sup>SM</sup>.

Nordic System for validation of alternative microbiological methods

**Nordval Validation (Ref. No.: 2003-20-5408-00019)** for all foods

Refer to NordVal validation for Petrifilm Staph Express Plate method details.

**AFNOR validated method in compliance with ISO 16140 in comparison to ISO 6888-1 (3M – 01/9 – 04/03A) and ISO 6888-2 (3M – 01/9 – 04/03B)**

**Scope of the validations:** All foods.

Use the following details when implementing the above Instructions for Use:

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1. FDA. 1998. Bacteriological Analytical Manual, 8th ed., (Revision A), Appendix 3.64
2. International Standards Organization, ISO 6887-1:1999, Microbiology – General guidance for the preparation of dilutions for microbiological examination
3. International Standards Organization, ISO 8261:2001, Milk and Milk Products – Preparation of test samples and dilutions for microbiological examination

AOAC is a registered trademark of AOAC INTERNATIONAL

Official Methods is a service mark of AOAC INTERNATIONAL







AFNOR is the trademark of Association française de normalization

NordVal is established as a Nordic System for validation of alternative microbiological methods (test kits).

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9 - 15 Chilvers Road  
Thornleigh, NSW 2120  
Australia  
1300 363 878

**3M****3M Microbiology**

St. Paul, MN 55144-1000  
Made in U.S.A. of globally sourced materials.

3M and Petrifilm are trademarks of 3M Company  
© 3M 2009. All rights reserved. 34-8704-4956-7

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3M Deutschland GmbH